



## **Offset Canes**

The **ProBasics Offset Cane** allows for easy height adjustment with a simple push of a button while the soft cushioned handle with convenient wrist strap improves comfort and grip. Available in three attractive finishes with a limited lifetime warranty.



Easy one-button height adjustment



## FEATURES

- Easy one-button height adjustment with a locking anti-rattle ring allowing for quick, safe adjustment
- Adjusts in 1-inch increments
- Comfortable, contoured handle
- Easy-to-clean aluminum frame
- Attractive powder-coated finish available in silver, black or bronze
- Bariatric model features reinforced tubing and a wider base for increased stability and support
- HCPCS\* Code: E0100

## SPECIFICATIONS

Offset: 29" – 38" Bariatric Offset: 30" – 39"
Offset: 7/8" Bariatric Offset: 7/8"
Offset: 3/4" Bariatric Offset: 3/4"
Offset: 300 lb Bariatric Offset: 500 lb
Offset: 12 oz Bariatric Offset:2 lb 5oz
Offset: 33.1" x 7.9" x 7.3" Bariatric Offset: 30.7" x 11" x 9.1"
Offset: 8.4 lb Bariatric Offset: 47.75 lb
Limited Lifetime on frame 90 days on all parts

ltem #	Description	Color	UPC	UOM
CNOFBK	ProBasics Offset Cane with Wrist Strap	Black	815067071135	10/cs
CNOFBZ	ProBasics Offset Cane with Wrist Strap	Bronze	815067071142	10/cs
CNOFSL	ProBasics Offset Cane with Wrist Strap	Silver	815067071159	10/cs
CNOFHD	ProBasics Bariatric Offset Cane with Wrist Strap	Black	815067071173	20/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.