



Transport Rollator

The **ProBasics Transport Rollator** combines the features of a rollator with a transport chair all in one unit. An individual can ambulate independently or be pushed safely by a caregiver. The footrests flip-down and the padded backrest is reversible depending on whether you are using it as a rollator or transport chair.

Equipped large 8" casters, easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch and a fold-strap for easy storage and transport. Designed to meet all your patients' fundamental needs and is height adjustable for a custom fit. Supports patient weights up to 250 pounds.



Flip-Down Foot Rests
and Water-Resistant Pouch



Easy-to-Operate
Ergonomic Hand Brakes



8-inch Wheels

FEATURES

- Lightweight aluminum frame with burgundy or blue finish
- Fold-down footrests and reversible padded backrest for easy conversion from rollator to transport mode
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy
Finish



Blue
Finish



SPECIFICATIONS

Overall Width:	24"
Width b/w Hand Grips:	19"
Seat Dimensions:	14.25" (W) x 13.25" (D)
Seat-to-Floor Height:	21.75"
Weight Capacity:	250 lb
Wheel Size:	8"
Product Weight:	22 lb
Carton Dimensions:	25" x 10.2" x 34.25"
Shipping Weight:	25 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLATBG	ProBasics Transport Rollator	Burgundy	815067071807	1/cs
RLATBL	ProBasics Transport Rollator	Blue	815067071814	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.