



## Two-Button Walker with Roll-Up Seat



**ProBasics Two-Button Folding Walker with Roll-Up Seat** provides patients with the convenience of a rollator and the functionality of a walker. Nylon seat hooks on frame for a secure seat. When not in use, it rolls up and stores on the frame allowing for a full stride. Features easy-to-use push button mechanisms allowing each side to operate independently for movement through narrow spaces. Supports patient weights up to 300 pounds.



Nylon seat hooks on frame for a secure seat



Roll seat up when not in use



Push button provides audible click when locked

### FEATURES

- Durable nylon roll-up seat that can be placed on handy hooks when not in use
- Adjustable height settings
- Easy-to-use push button mechanisms
- Front cross braces near top of walker to allow for full stride
- Each side operates independently to allow for easy movement through narrow spaces
- Vinyl contoured hand grip for added comfort
- Attractive, easy-to-clean aluminum frame
- Folds flat for transport and storage
- Weight Capacity: 300 lb.
- HCPCS\* Code: E0143 (with wheels)  
E0156 (seat)

### SPECIFICATIONS

Width Inside Hand Grips:	18"
Width Inside Base Legs:	22½"
Depth Opened at Base:	18½"
Depth Folded:	4"
Patient Height Range:	5'2" – 6'4"
Weight Capacity:	300 lb
Product Weight:	8.45 lb
Carton Dimensions:	21.7" x 6.3" x 31.5"
Shipping Weight:	10 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
WKA AW2BST	ProBasics Two-Button Walker with 5" Wheels and Roll-Up Seat	815067071098	2/cs

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.