


BED ACCESSORIES

Reduced Gap Full-Length Bed Rails #6034

#6034

- Welded steel construction with brown vein finish
- Self adjusting spring loaded crossbraces
- Four crossbars reduce entrapment hazards
- Overall length: 55"
- Overall height above deck: 13.5" raised - 4.5" lowered



Four Bar Half Rails #7035

- No gap style
- Welded steel construction with brown vein finish
- Fits all woven and spring style hospital beds
- No tools required
- Height 19", length 37"

Overbed Tables #3037

- Non Tilt
- Easily assembled
- Swivel ball bearing casters
- 30" Wood grain laminated top (30"w x 15"d)
- "H" base with brown powder coat finish
- Infinite height adjustment from 27" to 42"
- Supports up to 25 lbs. in lowest position

