

# Geo-Mattress with Wings<sup>®</sup>



Item #	Description
W7535-29	75"L x 35"W x 6"H center, 8"H side
W7835-29	78"L x 35"W x 6"H center, 8"H side
W8435-29	84"L x 35"W x 6"H center, 8"H side



**Stable, high-density angled bolster**, topped with proprietary Geo-Matt<sup>®</sup> segmented surface.



BEDS



SEATING



POSITIONERS



SKIN CARE

## Geo-Mattress with Wings<sup>®</sup>

*Safety – without compromise.*

Helps protect skin integrity while providing gentle prompting away from the edge of the bed.

- 5-year, non-prorated warranty
- Weight Limit: 400 lbs.

Flammability Protection: Certified conformance to NFPA 101 (Life Safety Code); Cal. Tech Bulletin #117; 16 CFR 1632 and 1633

### Primary Use:

- Prevention of Skin Breakdown


### Patient Indications:

- Patients at fall risk due to independent exits from bed
- General patient population
- Moderate to high risk for skin breakdown
- Decreased mobility
- Repositions frequently by self or caregiver


### Key Advantages:

- Edge-to-edge antishearing surface
- Extra-firm edge bolster
- Shearing protection at transfer opening
- Heel Slope<sup>®</sup>

### Packed with value-added features:




**Lowered, padded openings** make transfers easier and safer. Proprietary edge-to-edge Geo-Matt segmented surface provides shearing protection.



**Exclusive Heel Slope<sup>™</sup>.** Subtle taper redistributes load to pressure-tolerant lower legs. Delivers heel pressures 27% lower than flat designs.



**Fluid-resistant, wipe-clean cover** with ultra-low moisture vapor transmission rate (MVTR). Bacteriostatically-treated top, durable tri-laminate bottom. Hidden zipper design for easier, more thorough cleaning and superior infection control.



**Effective pressure redistribution** across the surface ensures low average of peak pressure at vulnerable sacrum and heels.

**SPAN**America  
Innovative Solutions.