

Be Independent™


SECURITY POLE

BECAUSE YOUR SECURITY POLE CAN GO ANYWHERE IN THE HOUSE



PRODUCT #1150





FEATURES & BENEFITS

- **Comfortable Cushion Grip** – Easy-to-sanitize closed cell foam provides non-slip surface.
- **Adjustable** – Fits ceiling heights 7'–10' ft.
- **Portable** – Main pole easily separates into two 53" pieces making it easy to transport.
- **Adaptable** – No need for wall mounts; can be installed in virtually any area of a room exactly where it is needed.

SPECIFICATIONS

Materials: Rust resistant zinc plated steel with durable powder coat finish

Fits Ceiling Heights: Range of 7'-10' ft.

Weight Capacity: 300 lbs

Weight of Product: 20 lbs

Package Dimensions: 54.00" L × 12.50" W × 3.25" H; 22 lbs

Case of 2 Dimensions: 54.00" L × 12.50" W × 6.5" H; 45 lbs

COMES IN TWO COLORS

 **BLACK #1150-B**

 **WHITE #1150-W**

WHAT IS IN THE BOX



**NO ADDITIONAL
TOOLS REQUIRED**



WEIGHT CAPACITY