

SmartCore®

Disposable Briefs

Ideal Micro-Climate for Skin Integrity Daytime/Nighttime

An Ideal Micro-Climate is achieved as the breathable sides provide adequate air-flow allowing heat to escape. Moisture is absorbed into the high capacity core keeping skin dry and comfortable.











Performance Advantage:

- **Longer Wear-time** – Clinical trials have proven that the brief can remain effective up to 6 to 8 hours without changing depending on fluid intake.
- **Fewer Changes** – brief holds multiple voids reducing the number of changes required and reducing physical strain on staff and caregivers.
- **Uninterrupted Sleep** – better sleep patterns through the night, positively impact overall health.

Product Features:

- **Peach Mat Construction** – a patented composite structure that guarantees superior performance in:
 - Skin Dryness
 - Neutralization of Urine pH
 - Odor Reduction
 - Inhibition of Bacterial Growth.
- **Breathable Sides** – increased fresh air flow eliminates heat and moisture build-up to promote skin dryness.
- **Kufguards®** – inner leg cuffs designed to channel fluids, contain high volume urge episodes and help contain bowel incontinence.
- **Soft, Cloth-like Outer Layer**
- **Micro-hook Closure Tabs** – grip the cloth-like outer layer to securely fasten the brief and allow multiple openings.
- **Gentle Leg Gathers** – for a secure fit.
- **Wetness Indicator** – dual, yellow lines turn blue when moisture is present, thus eliminating unnecessary changes and reducing costs.
- **Color-coded in 5 Sizes** – differentiation of size for simplifying caregiver duties and helping with proper size selection.
- **Latex-free**



Item #	Description	Size	Capacity Under Pressure*	HCPC	Packaging
2311	Small 	24 _ 32"	22.0 oz / 650 mL 	T4521	10-10's (100/cs)
2312	Medium 	32 - 44"	28.0 oz / 828 mL 	T4522	8-12's (96/cs)
2313	Large 	45 - 58"	34.0 oz / 1005 mL 	T4523	8-12's (96/cs)
2314	X-Large 	56 - 64"	34.0 oz / 1005 mL 	T4524	6-12's (72/cs)
2315	2XL 	60 - 80"	34.0 oz / 1005 mL 	T4524	4-8's (32/cs)

*C.U.P. is the amount of fluid absorbed by the product that cannot be forced out. Normal urine void is 8-12 ounces.

Free Samples and Consultation
Call Tranquility's Care Center
1-800-467-3224

© Tranquility is a registered trademark of Principle Business Enterprises, Inc.
© Copyright 2015 | 10TRQ-695-SL15

TRANQUILITY®
PREMIUM PROTECTION

www.tranquilityproducts.com
1-800-467-3224